

SteriTite® Sealed Containers

Universal Compatibility: SteriTite® is the ideal reusable packaging system for sterilization and transport of medical devices. It is durable, “virtually repair free”, more efficient and cost-effective than other containers on the market. SteriTite containers with MediTray® products are compatible with all current sterilization modalities meeting FDA, ISO and EU regulatory requirements. Anodized aluminum and stainless steel baskets are available to secure the contents.

Ideal for Standardization: SteriTite containers are sold as a complete unit with lid, base, and filter retention plates. The part numbers below represent our line of solid bottom containers for use in pre-vacuum steam and V-Pro low temp sterilizers. Perforated bottom units are available for gravity displacement steam and for STERRAD Sterilization; their part numbers are identified with the suffix "G". Request ID plaques for identification or barcodes for digital visualization.

Mini & Narrow

SC03M



SC02N



SC04NL



SC02M 12"x6"x2" (30x15x5 cm)

SC03M 12"x6"x3" (30x15x8 cm)

SC04M 12"x6"x4" (30x15x10 cm)

SC02N 18.5"x6"x2" (47x15x5 cm)

SC03N 18.5"x6"x3" (47x15x8 cm)

SC03NL 27"x6"x3" (68.5x15x10 cm)

SC04NL 27"x6"x4" (68.5x15x10 cm)

SC05NL 27"x6"x5" (68.5x15x12.7 cm)

Half Size



SC04H

SC04H 11"x11"x4" (28x28x10 cm)

SC05H 11"x11"x5" (28x28x13 cm)

SC06H 11"x11"x6" (28x28x15 cm)

SC08H 11"x11"x8" (28x28x20 cm)

Three Quarter Size



SC04Q

SC04Q 18.1"x11"x4" (46x28x10 cm)

SC05Q 18.1"x11"x5" (46x28x13 cm)

SC06Q 18.1"x11"x6" (46x28x15 cm)

SC08Q 18.1"x11"x8" (46x28x20 cm)

Full Size



SC04F

SC04F 23.5"x11"x4" (60x28x10 cm)

SC05F 23.5"x11"x5" (60x28x13 cm)

SC06F 23.5"x11"x6" (60x28x15 cm)

SC08F 23.5"x11"x8" (60x28x20 cm)

Long Size



SC04L

SC04L 27"x11"x4" (69x28x10 cm)

SC06L 27"x11"x6" (68x28x13 cm)

SC08L 27"x11"x8" (69x28x15 cm)

Wide Size



SC05W

SC05W 26.5"x16"x5" (67x41x13 cm)

SteriTite containers

Use disposable cellulosic filters for steam and polypropylene filters for low temp sterilization, as well as tamper evident seals and load cards with chemical indicator provided by Case Medical.



Customization

SteriTite® containers are the ideal choice for standardization

- **Universally compatible:** Cleared for use with all current sterilization modalities and medical devices that can be sterilized
- **Faster Throughput:** Validated for 5 to 8 minute dry times for terminal steam sterilization - reduce sterilization time and eliminate IUSS
- **Prevent strikethrough:** Our unique offset filtered vent pattern prevents cross contamination and ensures that the sterility of the contents is maintained
- **UDI ready:** Each container has a unique UDI barcode for seamless tracking, tracing, and direct access to the IFU
- **Ideal for loaner drop-in:** Our DIN-sized containers save time - just drop in your loaner or consignment tray to eliminate the need to wrap
- **Save money:** Lower cost per use than sterilization wrap and a lower cost of ownership than competitive containers


Basic customization with the MediTray® System

MediTray baskets and trays are designed to fit seamlessly in SteriTite containers in single layers or stacked for efficient use of container space. Select the right insert for sterilization, utilize our modular components, or consider one of our many customized options.




Customization with the MediTray® System


Case Medical offers modular, cost-effective, trays and components to customize and organize virtually any procedure tray. Our Web Sets are pre-assembled tray systems that provide the ultimate in instrument protection and fit seamlessly into our SteriTite container.




Cysto
YSC1003A2W



Da Vinci Xi
Endoscope Tray
YNB1009021H



Flexible Endoscope
YHA1024A9A



Amputation Tray
YRA1004D9C



Camera Case
YRA1253A12H



Cataract Set
YPA1002A3E